POWER ////



Easy-to-use controller with speed dial, battery indicator and hand controlled lever

GIVE YOUR MANUAL WHEELCHAIR A HEAVY DUTY BOOST

The Power Assist HD is designed to aid caregivers in pushing a manual wheelchair. The Power Assist HD is compatible for manual wheelchairs with widths of 16"-24" and can last up to 16km (10 miles) between charges.



20-24 Apollo Drive Hallam VIC 3803 PH: (03) 8770 9600 | FAX: (03) 9703 2720

www.pridemobility.com.au ACN 088 609 661 ABN 66 088 609 6611







SPECIFICATIONS

Maximum Speed ¹		6.4 km/h (4mph)
Cruising Range between Charges ¹	up to	10 miles (16km)
Maximum Gradient 6°		
Weight Capacity		182kg (400lbs)
Motor		24V 200W D.C.
Weight	• •	ight with battery) kg (battery pack)
Battery 2x	12V 12ah S	ealed Lead Acid
Charger		2A 24V output
Controller	PG	- S - DRIVE 45A
Fits Wheelchair Widths	40.6cm -	61cm (16" - 24")
Dimensions ²	Height: Length: Width:	38cm (15") 38cm (15") 25cm (9.8")

¹ Speed and range vary with user weight, terrain type, battery charge, battery condition and tyre pressure.

 2 Due to manufacturing tolerances and continued product improvement, this specification is subject to a variance of +/- 3%.



07 3265 4663 2/452 Bilsen Rd, Geebung 0LD 4034 sales@powermobility.com.au • powermobility.com.au



The information contained herein is correct at the time of publication; due to the manufacturer's commitment to constant improvement and development. They reserve the right to alter specifications without prior notice.



FEATURES

- Designed to work with manual wheelchairs with widths 16" - 24".
- Range of up to 16km (10 miles) between charges.
- Weight capacity of 182kg/400lbs (user and wheelchair combined).
- Maximum speed of up to 6.4km/h (4mph).
- Dual-wheel design for added traction and stability

